



PHOTOFRIN® (porfimer sodium) 75 mg for Injection

PHOTOFRIN® (porfimer sodium) is available exclusively through CuraScript, Inc.

OPTIGUIDE™ Fiber Optic (DCYL Cylindrical Diffuser Series)

The OPTIGUIDE™ DCYL200 series optical diffuser fiber is a light delivery system consisting of a 400µm core silica optical fiber, a proximal SMA-type connector, and a distal light diffusing tip. The fiber is designed to distribute light energy uniformly over the specified length of the diffuser tip.

The OPTIGUIDE™ DCYL200 series fiber is intended for single use only. All OPTIGUIDE™ DCYL200 series fibers are intended for use with DIOMED®, Coherent®, and Laserscope® PDT lasers.

CuraScript Item Number	Model Number	Diffuser Length/Description	Tip Outer Diameter	List Price (ea.)
245283	OPTIGUIDE™ DCYL210	1.0 cm Cylindrical	1.65 mm	\$895.00
245291	OPTIGUIDE™ DCYL215	1.5 cm Cylindrical	1.65 mm	\$895.00
245305	OPTIGUIDE™ DCYL220	2.0 cm Cylindrical	1.65 mm	\$895.00
245313	OPTIGUIDE™ DCYL225	2.5 cm Cylindrical	1.65 mm	\$895.00
245321	OPTIGUIDE™ DCYL250	5.0 cm Cylindrical	1.65 mm	\$895.00
255432	OPTIGUIDE™ DCYL710	1.0 cm Cylindrical	1.07 mm	\$895.00
255440	OPTIGUIDE™ DCYL725	2.5 cm Cylindrical	1.07 mm	\$895.00
255459	OPTIGUIDE™ DCYL750	5.0 cm Cylindrical	1.07 mm	\$895.00

All rights reserved.

PHOTOFRIN® (porfimer sodium) and OPTIGUIDE™ are registered trademarks of Concordia Laboratories Inc.

Pinnacle Biologics™ and the logo of Pinnacle Biologics™ are trademarks of Pinnacle Biologics, Inc.

PHOTOFRIN® is manufactured by Jubilant Hollisterstier General Partnership for Concordia Laboratories Inc. and distributed in the United States by Pinnacle Biologics Inc.

See more at: <http://www.photofrin.com/healthcare-professional-home/#sthash.FInUv0e7.dpuf>.

2015 Pinnacle Biologics Chicago, IL 60606

PINMAR002 1/2015 All rights reserved.

Ordering Information

To place your order for OPTIGUIDE™ fibers or for PHOTOFRIN®, please order through CuraScript, Inc.

Phone: (866) 268-5554

Customer Service: (877) 599-7748

Fax: (800) 862-6208

Website: www.curascriptonline.com

CuraScript, Inc. can be contacted:

Monday – Friday

8:30 a.m. – 7:00 p.m. Eastern Time

(Orders will be delivered the next business day)

**Announcing the
newest addition to
our fiber portfolio —
DCYL700 series is here!**