



PHOTOFRIN® (porfimer sodium) for Injection

PHOTOFRIN® (porfimer sodium) for Injection is available exclusively through CuraScript, Inc.

OPTIGUIDE™ Fiber Optic

The OPTIGUIDE® optical diffuser fiber is a light delivery system consisting of a 400µm core silica optical fiber, a proximal SMA-type connector, and a distal light diffusing tip. The fiber is designed to distribute light energy uniformly over the specified length of the diffuser tip.

The OPTIGUIDE® series fiber is intended for single use only. All OPTIGUIDE® series fibers are intended for use with 630 PDT lasers.

CuraScript SAP Material Number	Diffuser Length	Type	Tip Outer Diameter	List Price (Each)
10004619	1.0cm	Rigid	1.7mm	\$895.00
10004620	1.5cm	Rigid	1.7mm	\$895.00
10004621	2.0cm	Rigid	1.7mm	\$895.00
10004622	2.5cm	Rigid	1.7mm	\$895.00
10004623	5.0cm	Rigid	1.7mm	\$895.00
10004654	1.0cm	Flexible	1.2mm	\$895.00
10005192	2.5cm	Flexible	1.2mm	\$895.00
10004655	5.0cm	Flexible	1.2mm	\$895.00

PHOTOFRIN® and OPTIGUIDE® are registered trademarks of Concordia Laboratories Inc. Pinnacle Biologics™ and the logo of Pinnacle Biologics™ are trademarks of Pinnacle Biologics, Inc. PHOTOFRIN® is distributed in the United States by Pinnacle Biologics, Inc. US/PHO/NP/0018 (08/2024)

Pinnacle Biologics, Inc. Bannockburn, IL 60015 All rights reserved.

ORDERING INFORMATION

To order OPTIGUIDE® fibers or PHOTOFRIN®, please contact CuraScript, Inc.

Phone: (866) 268-5554

Customer Service: (877) 599-7748

Fax: (800) 862-6208

Website: www.curascriptonline.com

CuraScript, Inc. can be contacted:

Monday - Friday

8:30 a.m. – 7:00 p.m. Eastern Time

(Orders will be delivered the next business day)

Product Inquiries

Reports about safety, adverse events, side effects and product complaints can be made by calling 1-866-248-2039 or by emailing medicalinformation@advanzpharma.com

